

WATER BATH



- ▶ LED digital microprocessor controlled water bath with precision temperature control within +/- 0.5C (+/- 1F)
- ▶ 120V/60hz, 1440 Watt heating element provides quick heat-up and even heat (North American specs)
- ▶ Met Lab approved to UL standards 60601-1 and CSA #C22.2NO.601
- ▶ Exterior Dim: 84cm x 66cm x 20cm (33"x26"x8")
Interior Dim: 66cm x 56cm x 13cm (26"x22"x5")
- ▶ Side mounted drain with lockout for safe easy water exchange




WB-LRG-2

Digital Water bath

WB-LRG-2-EU

Digital Water bath
(220V/50hz)



Water Drain



Micro Processor Control

DRY THERMOPLASTIC OVEN



CO-SM-1

Small compact convection Thermoplastic oven fits most LT Thermoplastic masks with manual control knobs and timer

Exterior Dim: 12.5" x 20.75" x 15.5"

120V/60hz, 1750Watt, CE marked and NSF approved



CO-LRG-2

Larger convection Thermoplastic oven fits all LT Thermoplastic masks and ShoulderTech™. Thermoplastics digitally controlled

Exterior Dim: 24" x 28" x 17"

Interior Dim: 19.5" x 15.5" x 10"

120V/60hz, 1450Watt, UL, NSF and CSA approved

Ask about our other thermoplastic dry heat ovens to fit your specific requirements